

Catalog #	RK-071-09
Standard Size	125 RIA tubes
Radioactivity	This kit contains 1.5 µCi of Iodine 125
Range	1-128 pg/tube in 100µl of Standard Solution = 10-1280 pg/ml
Lowest Detection Limit	32 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Urotensin II (Rat)	100
	Urotensin II-Related Peptide	100
	Urotensin II (Mouse)	20
	Urotensin II (Human)	1
	Endothelin-1	
	(Human, Porcine, Rat Mouse, Canine, Bovine)	0
	CN-22 (Human, Porcine, Rat)	0
	CGRP (Rat)	0
	Adrenomedullin (1-50) (Rat)	0
	ANP (Rat, Rabbit, Mouse)	0
	Bradykinin	0
	Angiotensin II (Human)	0
	BNP-45 (Rat)	0
	BNP-32 (Rat)	0
	PAMP-20 (Rat)	0

Sequence pGlu - His - Gly - Thr - Ala - Pro - Glu - Cys - Phe - Trp - Lys - Tyr - Cys - Ile

Disulfide Bridge: Cys8 - Cys13