## Anti ACTH (1-23) (Mouse, Rat) Serum

Cat. No. Y352 Lot No. 190100513

**Description:** This antiserum was raised in a rabbit by immunization with a keyhole lympet hemocyanin (KLH) conjugate of synthetic ACTH (1-23) (mouse, rat) peptide. The product vial contains 50  $\mu$ L of the titled antiserum obtained by lyophilizing its 0.001M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohi- stochemistry or any other immunoreaction with the N-terminal portion of ACTH (mouse, rat).

Immunogen: Synthetic ACTH (1-23) (mouse, rat)-KLH conjugate Host: Rabbit

Amino Acid Sequence of ACTH (1-23) (mouse, rat) :

SYSMEHFRWG KPVGKKRRPV KVYPNVAENE SAEAFPLEF

**Product Form:** Lyophilized unpurified serum

Size:  $50 \mu L$ 

Reconstitution: Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN3 0.1%).

Storage: The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon reconstitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezingthawing should be avoided.

Suggested Working Dilution Range: 1:1,500 (final dilution ~1:10,500) for radioimmunoassay; 1:200-1,000 for immunohistochemistry (frozen or paraffin section). Optimal dilution should be determined by each laboratory for each application.

Specificity (based on enzyme imunoassay): ACTH (1-23) (mouse, rat) 100%, ACTH (mouse, rat) 100%

Positive Control (immunohistochemistry): Rat pituitary gland

Species Tested: Rat

## REFERENCES:

1) J. Drouin, M. Chamberland et al., Structure of the rat pro-opiomelanocortin (POMC) gene. FEBS Letters. 193:54-58, 1985

## FOR RESEARCH USE ONLY

< Distributed by >

SCETI K.K.

DF Kasumigaseki Place, 3-6-7 Kasumigaseki, Chiyoda-ku, Tokyo 100-0013 Japan

URL: http://www.sceti.co.jp/export/e-mail: exp-pet@sceti.co.jp

<Manufacturer>

Yanaihara Institute Inc.

2480-1 Awakura, Fujinomiya-shi, Shizuoka 418-0011 JAPAN