

Anti-HIV-1 p55 antibody

Immuned Animal: Rabbit

Polyclonal antiserum

65-013 50 µl, 65-014 250 µl

HIV-1 Gag p55 is a precursor protein of several proteins that form the core structure of AIDS virus, indispensable to its reproduction. P55 is digested by HIV-1 protease, first into intermediate products p41 and p15. Then p41 is digested into matrix protein p17 and capsid protein p24. Protein p15 is further digested into nucleocapsid protein p7 and to p6 and p1, both of unknown function (1).

The product is prepared by immunizing rabbit with recombinant protein p55 which was over-expressed in *E. coli* with a plasmid carrying the Gag p55 coding region of HIV-1 virus, and highly purified by several steps of chromatography (2).

Using this antiserum in Western blotting, in addition to the band of 55kD for HIV-1 p55, bands of 15 kD, 17 kD and 24 kD were observed in the extract of the AIDS virus infected cells (Fig. 1).

Usage

1) It can be used in Western blotting or ELISA for HIV-1 p55 detection and titration.

Specifiication

Form: 0.09% sodium azide added to the antiserum.

Storage: 4°C

Reference:


1. Freed EO, Virology 251:1-15 (1998) Review
2. Saito A, et al., Microbiol Immunol 39:473-483 (1995)

Fig. 1 Detection of HIV-1 p55 by Western blotting using the p55 antibody.

Lane 1, Extract of MT4 cells.

Lane 2, Extract of MT4 cells infected with HIV-1 (LAI strain).

The antiserum was diluted 2000 fold before use.



<Distributed by >: SCETI K.K. DF Kasumigaseki Place, 3-6-7 Kasumigaseki,
Chiyoda-ku Tokyo 100-0013 JAPAN
Tel: +81-3-5510-2347 Fax: +81-3-5510-0134
E-mail: exp-pet@sceti.co.jp URL: www.sceti.co.jp/export/

<Manufactured by> : BioAcademia, Inc. 7-7-18 Saito-Asagi, Ibaraki, Osaka 567-0085, JAPAN